

# MDI Study Viewer

The MDI Study Viewer allows you to use the functionality of our digital image management system (AGFA PACS) to electronically view films in your surgery.

You will benefit by being able to use an online diagnostic imaging application, whereby patient images, past studies and reports are securely available, as quickly as your Internet connection allows.

Functions include viewing tools available to the radiologist, as well as annotation and exportation facilities. Direct transfer of patient data to CD/DVD-R is provided, allowing patients to take a copy of their films/reports to other doctors, if required.

MDI stores an electronic record of images (normally for seven years) and these are available to you for retrieval.

The MDI study viewer can usually be added to your existing IT system. Minimum system requirements are:

- Windows 2000, XP or VISTA
- Pentium 4 - 1.5 GHz or equivalent
- 512Mb of RAM
- 200Mb of available hard disk space
- VGA compatible video adapter and monitor
- CD-ROM or DVD drive
- Keyboard and Mouse or equivalent pointing device

We hope you will take advantage of this useful tool for radiology image viewing and management. Impax application is fully complimentary.

Please note: Patients will still receive hard copy films unless otherwise advised.